

DR 200™

No Rinse Disinfectant Sanitizer



Description

DR 200™ is a concentrated disinfectant and sanitizer that requires no rinsing. It is designed specifically for hospitals, food processing plants, dairies, restaurants, bars, animal quarters, kennels and any institution where disinfection and sanitization are of prime importance. Depending on its dilution, **DR 200™** can either disinfect or sanitize pre-cleaned inanimate hard surfaces such as walls, floors, sink tops, tables, chairs and bed frames. **FOR PROFESSIONAL USE ONLY.**




Benefits

- Controls mildew and bad smells.
- pH-neutral.
- Effective against several types of viruses.

pH 7

Directions

1. Wet all surfaces thoroughly with the solution using a cloth, mop, microfiber towel, spray bottle or spray gun. **2.** Allow about 10 minutes before removing excess fluid. **3.** Let air dry. Always clean before disinfecting hard, inanimate surfaces. Prepare a fresh solution for each use.



When using a spray bottle, adjust the tip so that the disinfectant is sprayed in droplets and not in mist.

- Consult an InnuScience representative for more information on the dilutions, dosages or equipment to be used according to your needs.
- Consult a representative about our specialized training.

Dilution Guideline

Disinfection - Let stand for at least 10 minutes

	600 ml	0.55 %	3.3 ml / 600 ml	1 : 180
	4 L	0.55 %	22 ml / 4 L	1 : 180

Mildewstat - Let stand for at least 10 minutes

	600 ml	0.81 %	4.9 ml / 600 ml	1 : 123
	4 L	0.81 %	32.5 ml / 4 L	1 : 123

Sanitization - Let stand for at least 60 secondes

	600 ml	0.2 %	1.2 ml / 600 ml	1 : 500
	4 L	0.2 %	8 ml / 4 L	1 : 500

Products Codes and Packagings

15666

4 L

2 Units